



#### Personal laser device

No external power cords or laser firing pedals are required during use, maximizing mobility and portability. (You can also connect a laser launch pedal, there are two modes to choose freely) The built-in long-life lithium battery can last for one month on a single charge.

# Dual wavelength: 650nm+980nm Easy to operate, every one can learn to use it.

Start up without waiting, use it immediately;

More than 10 built-in clinical indication parameter templates can automatically save your settings, which is simple and saves your time.

## Pay attention to every little detail of use -

---- The handle is removable and also autoclavable

--- Handpiece has power adjustment button which makes it easy to operate

---- Laser switch buttons ensure maximum safety

### Fiber optic tip models of Dental diode laser

Non-initiated	Models	ф	Length	Intended
141	F4-4	400um	4	Surgery
171	F4-7	400um	7	Periodontz
1 9	F4-9	400um	9	Periodonta
1 14 h	F2-14	200um	14	Periodonta
1 20	F2-20	200um	20	Root canna
mm				

Pre-initiated	Models	Ф	Length	Intended
141	PF4-4	400um	4	Surgery
1 7 h	PF4-7	400um	7	Periodontal
1 9 h	PF4-9	400um	9	Periodontal

optical fiber tips are autoclavable consumables

# **Optional**



Whitening-tip

Teeth Whitening



Physiotherapy-tip

Wound healing
Burning Mouth
Dentin Hypersensitivity